
PRODUCT SHEET



Eppendorf ThermoMixer® F2.0

The Eppendorf ThermoMixer F2.0 heats and shakes up to 24 x 2.0 ml tubes.

We also offer Eppendorf ThermoMixers® F2.0 for a wide range of other tube sizes like microplates, 0,5ml tubes or 1,5ml Tubes. Please inquire!

TECHNICAL DATA:

Capacity: 24 x 2,0 ml Tubes*
Speed: 300 - 1.500 rpm
Temperature: 4°C above RT to 100°C
Dimensions (WxDxH): 206x304x170 mm
Weight: 6,3 kg

* Note: We also offer ThermoMixers® F2.0 for all other tube sizes.

SERVICE:



Calibration

SERVICE:



reCAL Exchange
Program