

---

## PRODUCT SHEET



### Thermo Heratherm IMH180 Incubator

The Heratherm IMH180 Incubator features a volume of 178 liters and a temperature range of up to 105 °C. The adjustable fan speed ensures a high temperature uniformity throughout the chamber.

We offer a wide range of incubators in different sizes. Please get in touch to discuss your requirements.

### TECHNICAL DATA:

---

Volume: 178 l

Forced air circulation: Yes

Temperature: 5 °C above RT to 105 °C

Dimensions exterior (WxDxH): 640 x 920 x 738 mm

Dimensions interior (WxDxH): 464 x 408 x 543 mm

Weight: 70 kg

### SERVICE:



Calibration