
PRODUCT SHEET



Thermo Heracell 150i CO2-Incubator

The Thermo Heracell 150i is a high-quality CO₂ incubator with a precise CO₂ and temperature regulation and fast recovery time.

TECHNICAL DATA:

Volume: 150 l

CO₂ range: 0% - 20%

Temperature range: 3 °C above RT to 55 °C

Humidity: <90% at 37 °C

Dimensions exterior (WxDxH): 637 x 782 x 867 mm

Dimensions interior (WxDxH): 470 x 530 x 607 mm

Weight: 70 kg

SERVICE:



Calibration

ANCILLARIES:



Temperatur Datalogger
EBI 310