
PRODUCT SHEET



Liebherr LGUex 1500 Freezer (-20°C)

The Liebherr LGUex 1500 is our reliable, compact laboratory freezer.

TECHNICAL DATA:

Volume: 139 l

Temperature: -9°C to -26°C

Cooling system. Static

Fluctuation: max. 3,3°C

Gradient: max. 6,6 °C

Dimensions exterior (WxDxH): 600 x 615 x 820 mm

Weight: 42 kg

SERVICE:



Calibration

SERVICE:



Temperature Mapping

SERVICE:



IQ/OQ

ANCILLARIES:



Temperatur Datalogger
EBI 310 and
EBI310-TPX220