

---

## PRODUCT SHEET



### Liebherr MKv 3910 Pharmacy Refrigerator (+5 °C) (361 liters)

The Liebherr MKv-3910 ist a high-quality pharmacy refrigerator ideal for the use in clinical trials. It features a large 361 liter storage volume and a fixed temperature setting at 5 °C which eliminates any accidental manipulation. An optical and acustical alarm is given, when the temperature breaches the temperature limits +2 °C or +8 °C.

#### TECHNICAL DATA:

---

Volume: 361 l  
Temperature: +5 °C (fixed setting)  
Cooling system: forced air  
Fluctuation: max. 3,6 °C  
Gradient: max. 4,1 °C  
Dimensions exterior (WxDxH): 601 x 618 x 1840 mm  
Weight: 67 kg  
Pharmaceutical storage: meets DIN 58345

#### SERVICE:



Calibration

#### SERVICE:



Temperature Mapping

#### SERVICE:



IQ/OQ

#### ANCILLARIES:



Temperatur  
Datalogger EBI 310