

LIEBHERR

MKv 3910

Gross capacity, total	361 l
Product dimensions (H/W/D)	184 / 60.1 / 61.8 cm
Energy consumption in 24h	0.846 kWh / 24h
SmartMonitoring	prepared

> Get help for service



General product information

Product group output channels	Pharmacy refrigerator with Comfort controller
GTIN	9005382238233

Output and consumption

Energy consumption in 24h	0.846 kWh / 24h
Noise level	48 dB

Net capacity, total	280 l
Gross capacity, total	361 l
Refrigerant	R 600a
Voltage	220-240 V ~
Frequency	50 Hz
Connection rating	1.5 A

Dimensions and weight

Product dimensions (H/W/D)	184 / 60.1 / 61.8 cm
Insulation	40 - 40 mm
Weight (with packaging)	73.00 kg
Weight (without packaging)	67.00 kg

Control and functions

Control	Electronic control system
Control unit	Buttons
Malfunction: Warning signal	optical and acoustical
potential-free contact	✓
Power failure alarm	immediately upon power failure for min. 12h
Maximum fluctuation	3.6 °C
Gradient	4.1 °C
SmartMonitoring	prepared

Design and features

Side wall material	steel
Colour	white
Door/Cover material	steel
Lock	fitted
Self-closing door	✓

Door hinges	right interchangeable
Material of interior containers	Plastic white
Type of castors	—
Door design	HardLine
Handle	Ergonomic slimline handle with opening mechanism

Please note: despite careful data management, we reserve the right to make technical changes; and discrepancies may occur between the image or text content and the original appliance.

LIEBHERR