
PRODUCT SHEET



Liebherr MKUv 1610 Pharmacy Refrigerator (+5 °C) (141 liters)

The Liebherr MKUv 1610 is a compact high-quality pharmacy refrigerator ideal for the use in clinical trials. It features a 141 liter storage volume and a fixed temperature setting at 5 °C which eliminates any accidental manipulation. An optical and acustical alarm is given, when the temperature breaches the temperature limits +2 °C or +8 °C.

TECHNICAL DATA:

Volume: 141 l
Temperature: +5 °C (fixed setting)
Cooling system: forced air
Fluctation: max. 1,3 °C
Gradient: max. 2,8 °C
Dimensions exterior (WxDxH): 601 x 618 x 820 mm
Weight: 42 kg
Pharmaceutical storage: meets DIN 58345

SERVICE:



Calibration

SERVICE:



Temperature Mapping

SERVICE:



IQ/OQ

ANCILLARIES:



Temperatur Datalogger
EBI 310