

---

## PRODUCT SHEET



EBI310-TPX220 Temperature-Datalogger for -80 °C

The EBI310-TPX220 temperature-datalogger features an exterior sensor and is ideal for temperature monitoring at ultralow temperatures (-70 °C/-80 °C).

### TECHNICAL DATA:

---

Temperature range: -200 °C - +250 °C

Precision: ±1.0 °C (-100 °C...-20 °C) / ±0.2 °C (-20 °C...+60 °C)

Data memory: 120.000 readings

#### Features:

- Tamper proof
- Easy to use through automatic PDF report generation, direct USB readout and configuration
- Minimize training costs and user errors through simple, easy to use design
- Eliminate costly and time consuming IQ/OQ procedures through standard reporting process
- Reduce IT costs due to the standardized software
- GMP, pharmaceutical law compliant; EN 12830, ATP, VO (EG) 37/205 compliant 21 CFR Part 11 conform (Winlog.basic)

We offer several services for our loggers:

- Calibration of the logger and Programming Service: Logger is pre-set according to customer's specifications
- Logger Exchange Program: Annual exchange of the old logger with a freshly calibrated logger.

SERVICE:



Calibration

SERVICE:



reCAL Exchange  
Program

SERVICE:

10100  
00101  
10100

Logger Programming  
Service