
PRODUCT SHEET



Compact 5 Patientenmonitor

Kompakter Multiparameter Patientenmonitor mit großem 7" TFT Farbbildschirm.

TECHNISCHE DATEN:

Misst NIBP, EKG, Pulsfrequenz, Respiration, SpO₂, Temperatur (optional)
 ST - Streckenanalyse
 EKG - Kurve 1 - kanlig
 Alarmfunktion aller Parameter
 Maße (BxTxH): 238 x 155 x 175
 Gewicht: 4,5 kg

SERVICE:



Annual Safety
Check

SERVICE:



Calibration

SERVICE:



reCAL Exchange
Program

SERVICE:



User Training

ZUBEHÖR:



Klebeelektroden