
PRODUCT SHEET



Ergoline ergoselect 100 P Ergometer

The classic ergometer for exercise ECG tests. Free programming of exercise protocols, intuitive operation and universal interfaces (digital, analog, remote start) for connecting to all commercially available PC-based ECG systems and ECG recorders set this ergometer apart. Multiple options and a comprehensive range of accessories are available.

TECHNICAL DATA:

Load: 6-999 Watt, speed-independent
Speed range: 30-130 rpm
Brake system: microprocessor controlled eddy current brake
Handlebar adjustment: 360° inclination
Saddle height adjustment: continuous mechanic
Patient Body height: 120 - 210 cm
Patient weight (max): 160 kg
Accuracy: according to DIN VDE 0750-238
Dimensions (LxWxH): ca. 90 × 46 × 133 cm
Weight: 54 kg

SERVICE:



Annual Safety Check

SERVICE:



Calibration

SERVICE:



User Training