

---

## PRODUCT SHEET



### Omron M500 Blood Pressure Monitor

The M500 Intelli IT is a clinically validated, Bluetooth enabled, blood pressure monitor.

Take the hassle out of cuff placement on your upper arm thanks to Omron's unique Intelli Wrap Cuff - giving accurate blood pressure measurements in any position around the upper arm.

The Omron M500 Intelli IT can take three consecutive blood pressure measurements at 30 second intervals, then displays the average to give you a realistic indication of blood pressure over time.

### TECHNICAL DATA:

---

Device Type: Upper Arm Blood Pressure Monitor

Cuff Type: Intelli Wrap Cuff 22-42 cm

AFib Detection: Yes

Memory: 2 Users x 100 Readings

Cuff Wrap Guide: Yes

Averaging Function: Yes

Irregular Heartbeat Detection: Yes

Body Movement Detection: Yes

Validation: Clinical Validation, Diabetic Validation, Pregnancy Validation

Intellisense Technology: Yes

Dimensions (WxDxH): 191 x 120 x 85 mm

Weight: 460 g

### SERVICE:



Calibration

### SERVICE:



reCAL Exchange  
Program